

MILK FATTY ACID ANALYSIS

FOSS

DR. DANIEL SCHWARZ, FOSS, DENMARK AND OLE RASMUSSEN, FOSS NORTH AMERICA

09 September 2019, Dallas, Texas, USA



1. INTRODUCTION

EVOLUTION OF UNDERSTANDING OF MILK FATTY ACIDS

FOSS

1897

- Question of the origins of milk fat (wholly from diet or synthesised by animal)? (Jordan and Jenter, 1897)

1947

- Leading theory: Short-chain fatty acids arise from degradation of oleic acid (Hilditch, 1947)

1951

- De novo synthesis of short-chain fatty acids proven (Popjak et al., 1951)

1990

- Detailed determination of origin of fatty acids completed (e.g., Palmquist, 2006)

2000s

- Fatty acid profiling (chain length, degree of saturation, major fatty acids) using FTIR technology (FOSS) → Practical applications: Focus on milk processing/dairy product aspects

2010s

- Practical applications: Fatty acid profiling with focus on dairy cow nutrition (Visiolait project)

FATTY ACID CALIBRATIONS AVAILABLE SINCE 2007

FOSS


ANALYTICS BEYOND MEASURE

- Degree of unsaturation
 - Saturated Fatty Acids (SFA)
 - Mono Unsaturated Fatty Acids (MUFA)
 - Poly Unsaturated Fatty Acids (PUFA)
 - Trans FA

- Chain length
 - Short Chain Fatty Acids (SCFA): $C_{4:0}$, $C_{6:0}$, $C_{8:0}$, $C_{10:0}$
 - Medium Chain Fatty Acids (MCFA): $C_{12:0}$, $C_{14:0}$, $C_{16:0}$
 - Long Chain Fatty Acids (LCFA): $C_{18:0}$, $C_{18:1}$, $C_{18:2}$

Major fatty acids

- $C_{14:0}$
- $C_{16:0}$
- $C_{18:0}$
- $C_{18:1}$




Package I

Package II


FREE FATTY ACIDS VS FATTY ACIDS

FOSS

Structural model of C12:0




Molecular structure of triglyceride



2. FATTY ACID UNITS

FAT AND FATTY ACIDS



→ Natural correlation between fat and fatty acids


| Fat (%) | Fatty acid (group), milk basis | Fatty acid (group), fat basis |
|---------|--------------------------------|-------------------------------|
| 1 | 0.1 | 10.52 |
| 2 | 0.2 | 10.52 |
| 3 | 0.3 | 10.52 |
| 4 | 0.4 | 10.52 |

according to Eskildsen et al., 2014


FATTY ACIDS – UNITS MILK BASIS VS FAT BASIS

FOSS

Milk basis




Fat basis



3. FATTY ACID ORIGIN PACKAGE

FATTY ACID ORIGIN GROUPING

FOSS



- Calibrations/prediction models based on natural material (i.e. raw milk) only
- Development of global models
 - Samples from around the globe included (reference and spectra data) → robustness
 - Variation in cow breeds and different seasons covered → robustness
- Regular updates of models (as industry evolves)
- Variation of results in samples more important than number of samples




VALIDATION OF DEVELOPED MODELS

FOSS

- Calibration development done on separate sample set; no calibration samples included!
- 303 samples available in total: 84 individual cow milk samples, 219 bulk tank/pen samples
- Geographical areas: 3 sites in USA, 2 in Canada, one each in Denmark, Germany, Japan, the Netherlands, Norway, Spain, Sweden, and the UK

VALIDATION – RESULTS ON BULK TANK/PEN SAMPLES

FOSS






ANALYTICS BEYOND MEASURE

VALIDATION – RESULTS ON INDIVIDUAL COW MILK SAMPLES

FOSS

A



4. QUALITY ASSURANCE

QUALITY ASSURANCE PROCEDURE

FOSS

1. Regular standardisation of MilkoScan
2. Adjustment against reference method or based on calibration samples
3. Monitoring of instrument based on check samples



FATTY ACIDS – CALIBRATION SAMPLES

FOSS

Unit: g FA/100 g milk

| Fatty acid origin groups | | |
|--------------------------|-------|-----------|
| De novo | Mixed | Preformed |
| 1.39 | 0.69 | 1.04 |
| 1.82 | 0.90 | 1.36 |
| 2.21 | 1.09 | 1.65 |
| 2.61 | 1.29 | |
| 2.97 | | |

Milk

**No true FA variation –
Not fit for purpose**

$$\frac{\text{Saturic Fatty Acid}}{100g \text{ Total Fatty Acids}}$$

Fat basis

| % fat |
|-------|
| 3.12 |
| 4.07 |
| 4.95 |
| 5.85 |
| 6.66 |

Unit: g FA/100 g TFA

| Fatty acid origin groups | | |
|--------------------------|-------|-----------|
| De novo | Mixed | Preformed |
| 46.9 | 23.2 | 35.2 |
| 46.9 | 23.2 | 35.2 |
| 46.9 | 23.2 | 35.2 |
| 46.9 | 23.2 | 35.2 |
| 46.9 | 23.2 | 35.2 |


EXAMPLE

DEVELOPMENT OF CALIBRATION SAMPLES

FOSS

- Bulk tank milk from different farms
- 14 samples offering true variation in FA composition

→ Excellent correlation between FTIR and reference results



FATTY ACID ORIGIN – DATA COLLECTION

FOSS


- Each month 14 calibration samples with GC results were prepared
 - Natural milk from mother nature (Adulteration not recommended)
 - Samples selected to have as wide a matrix variation as possible
 - Minimal covariance between different fatty acids
- Sample sets send to multiple instrument
 - MSC 7 (3 instruments – 2 different sites)
 - FT plus (3 instruments – Same site)
- First sample set tested end of March 2019 and last sample set tested early August (5 sets in total)
 - Time period between first and last adjustment is 4.5 months

TYPE OF ADJUSTMENTS – OPTIMAL ACCURACY

FOSS

- 70 samples collected and scanned over 4.5 months period
- 5 sets of 14 calibration samples
- **Experiment:**
 - Adjustments of Slope/Intercept every set (Months)
 - One adjustment used over entire time period


TYPE OF ADJUSTMENTS



TYPE OF ADJUSTMENTS

De novo FA (GC) VS Instruments

◆ De novo FA---MSC_7_A---S/I both/AllSets ■ De novo FA---MSC_7_A---S/I both/BySet



- Conclusions:
 - A single adjustment can be used for a long period of time
 - No need to adjust the calibration monthly
 - Developing guidelines for how best to manage calibration adjustments


TRANSFERABILITY – INSTRUMENT VS INSTRUMENT

FOSS

- The same 5 sets of calibration samples have been tested on 6 instruments
- Instruments known as **MSC_7_A** is considered the master instrument
- The master instrument is calibrated against reference chemistry (GC lab data)
- The other instruments are calibrated to the master instrument
- The room mean standard deviation of the differences has been calculated as follows

De novo FA---MSC_7_A VS Instruments


- ◆ De novo FA---MSC_7_C---S/I---BySampleSet
- De novo FA---MSC_7_B---S/I---BySampleSet
- ▲ De novo FA---FT_Plus_A---S/I---BySampleSet
- × De novo FA---FT_Plus_B---S/I---BySampleSet
- ✱ De novo FA---FT_Plus_C---S/I---BySampleSet



TRANSFERABILITY – INSTRUMENT VS INSTRUMENT

Residuals --- De novo FA

- ◆ De novo FA---MSC_7_C---S/I---BySampleSet
- De novo FA---MSC_7_B---S/I---BySampleSet
- ▲ De novo FA---FT_Plus_A---S/I---BySampleSet
- × De novo FA---FT_Plus_B---S/I---BySampleSet
- ✱ De novo FA---FT_Plus_C---S/I---BySampleSet



- The following Root Mean Standard deviation of differences has been calculated when comparing other instruments to the master.

| | | |
|---|-----------|-------|
| • | MSC_7_B | 0.007 |
| • | MSC_7_C | 0.009 |
| • | FT_Plus_A | 0.008 |
| • | FT_Plus_B | 0.011 |
| • | FT_Plus_C | 0.029 |

CONCLUSIONS ON TRANSFERABILITY – INSTRUMENT VS INSTRUMENT

FOSS


- The instrument to instrument agreement is significantly better than comparing to reference data (GC lab data)
- Satellite instruments can be calibrated to a master instrument.
- Master instrument must be calibrated against reference data (GC lab data)

5. POTENTIAL APPLICATIONS

HERD LEVEL



FAT AND FATTY ACID PROFILE

Fat naturally produced or through fat supplements?



REAL LIFE EXAMPLE FAT AND PROTEIN VS DE NOVO FA

FOSS




→ 200 dairy farms,
1 test per farm

High **fat** % associated with high de novo contents
→ increased function of rumen as well as production of volatile fatty acids

High **protein** % associated with high de novo contents
→ Increased microbial fermentation as well as microbial protein synthesis

→ Opportunity for dairy farmers to increase revenue and profit

REAL LIFE EXAMPLE: HERD MANAGEMENT



What happened?

- Cows started mobilizing
- De novo synthesis went down

Reason?

Change in silage quality
(more fibre/less digestible)

Idea:

Changes in fatty acid profile can be noticed a few days before milk or fat yield start to decrease

→ Dairy farmer can react EARLIER and save \$\$\$

COW LEVEL

INDIVIDUAL COW FA DATA – EXAMPLE

FOSS

Pen or group A

| FA level | Fatty acid group | | |
|----------|------------------|----------|-----------|
| | De novo | Mixed | Preformed |
| Low | 7 (35%) | 1 (5%) | 0 (0%) |
| Medium | 12 (60%) | 18 (90%) | 8 (40%) |
| High | 1 (5%) | 1 (5%) | 12 (60%) |

Pen or group B

| FA level | Fatty acid group | | |
|----------|------------------|----------|-----------|
| | De novo | Mixed | Preformed |
| Low | 6 (15%) | 2 (5%) | 23 (59%) |
| Medium | 14 (36%) | 26 (67%) | 9 (23%) |
| High | 19 (49%) | 11 (28%) | 7 (18%) |

DOCUMENTATION

FOSS

A WHITE PAPER FROM FOSS:

Fatty acid profiling according to origin for optimising feeding and management of dairy cows – a new approach

By: Dr. Daniel Schwarz
fossanalytic.com
November 2018, BB

FOSS

A WHITE PAPER FROM FOSS:

The new FOSS fatty acid origin package – Basics behind the prediction models

By: Dr. Daniel Schwarz, Maria Rosenberg Bak and Dr. Per Waaben Hansen, FOSS Denmark
fossanalytic.com • November 2018, BB


FOSS

Application Note

AN 5465

Rev. 2

MilkoScan™ 7 RM/ FT+ / FT 6000
Fatty Acid Origin Package



It is well-known that changes in both milk fatty acid composition and total fat content occur depending on the cow's diet. This Application Note describes three prediction models of fatty acid groups, which are based on the actual origin of fatty acids: de novo, mixed and performed.

Application Note 5465 / Rev. 2

FOSS

Fatty Acid Origin

A package of prediction models to help dairy farmers optimize dairy-herd feeding.




Milk fatty acid composition is clearly associated with the diet of dairy cows which, in turn, opens up the opportunity to improve herd efficiency via more effective feeding strategies. Specifically, the milk fatty acid origin profile of bulk tank samples allows effective evaluation of the feeding status of the entire lactating herd for improved fat and protein content in milk deliveries.

The Fatty Acid Origin package allows you to determine the fraction between the fatty acids according to their origin as guideline for feeding programs and management of cows. It consists of three prediction models: 1) de novo fatty acids, 2) mixed fatty acids, and 3) preformed fatty acids. The specific fatty acids included in each prediction model are shown in Fig. 1. overleaf.

The new models were developed based on full spectra data and gas chromatography reference data from multiple countries on three continents (Europe, North America, Asia). All data used originated from raw milk samples only (both bulk tank and individual cow milk samples) to reflect routine conditions in milk-testing laboratories.

Fat and fatty acids

- Fat in milk mainly appears in form of a complex combination of lipids called triglycerides
- Triglycerides consist of esters of three different fatty acids with one molecule of glycerol (see molecular structure)
- Fat content as well as the fatty acid composition are associated with animal-related factors and the management and feeding of the animal



ANALYTICS BEYOND MEASURE

ANALYTICS BEYOND MEASURE

A MESSAGE TO TAKE HOME

FOSS



Accurate, reliable, transferable, and stable generation of milk fatty acid profiles



New possibilities to improve feeding management



The power of co-creation



das@foss.dk



@SchwarzD123
@FOSSAnalytical



www.linkedin.com/in/daniel-schwarz84
www.linkedin.com/company/6750/